



07/05/2021

ASPIRE REHAB RS/RX TIS WHEELCHAIR RANGE CONSOLIDATION!

The Aspire Rehab RS TIS Wheelchairs will be discontinued, and we will continue only with Rehab RX models. Aside from colour, the technical differences between the RS & RX models are purely the seat and back cushioning... see below for details.

This consolidation enables us to hold more stock of the now fewer models (RX).

LIST and wholesale pricing of the RX will be adjusted slightly down to suit.

There is no change to the Aspire Rehab RX JNR, as it has never had an RS model variant.

KEY DIFFERENCES

- o **Backrest** – The RX has a multilayer foam contoured backrest with built up postural support and moderate lateral supports, while the RS has a standard foam, mildly shaped backrest cushion.
- o **Seat** – The RX has a multilayer foam seat cushion with integrated positioning skeleton as well as adductor and abductor contouring for enhanced sitting posture, while the RS has a standard foam seat cushion with limited shape.



Aspire Rehab RS

MWS449600 – Rehab **RS** 410mm – **Discontinued!**
 MWS449610 – Rehab **RS** 460mm – **Discontinued!**
 MWS449620 – Rehab **RS** 510mm – **Discontinued!**



Aspire Rehab RX

MWS449710 – Rehab **RX** 410mm – **Continues!**
 MWS449720 – Rehab **RX** 460mm – **Continues!**
 MWS449740 – Rehab **RX** 510mm – **Continues!**
 MWS449630 – Rehab **RX JNR** 320mm – **Continues!**

