

DVA PODIATRY PRODUCT LIST

Quality Equipment Enterprise Solutions Home & Community







DVA Podiatry Product List



Aidacare Locations

VICTORIA

MELBOURNE NORTH
82 Albion St,
Brunswick East VIC 3057
Tel: (03) 9981 2100
Fax: (03) 9384 2088

MELBOURNE SOUTH
313 Warrigal Rd,
Cheltenham VIC 3192
Tel: (03) 9584 1644
Fax: (03) 9584 0562

GEELONG
74-76 Fyans St,
Geelong VIC 3220
Tel: (03) 5223 2263
Fax: (03) 5223 2269

BENDIGO51 Bridge St,
Bendigo VIC 3550
Tel: (03) 5441 8213
Fax: (03) 5441 7762

MILDURA 433 San Mateo Ave, Mildura VIC 3500 Tel: 0400 938 062 Fax: (03) 5021 1845

NEW SOUTH WALES

\$YDNEY SOUTH
4 Noonan PI,
Ingleburn NSW 2565
Tel: (02) 8706 2300
Fax: (02) 9618 5111

SYDNEY NORTH
Unit 2/4 Stanton Rd,
Seven Hills NSW 2147
Tel: (02) 8706 2330
Fax: (02) 9624 4256

347 Urana Rd, Lavington NSW 2641 Tel: (02) 6040 3444 Fax: (02) 6040 3222

ORANGE
Unit 1/12-16 Ash St,
Orange NSW 2800
Tel: 0437 134 295

TAMWORTH
Unit 3/248 Peel St,
Tamworth NSW 2340
Tel: (02) 6761 2018
Fax: (02) 6761 2090

PORT MACQUARIE
31 Jindalee Rd,
Port Macquarie NSW 2444
Tel: (02) 6581 0018
Fax: (02) 6581 0038

COFFS HARBOUR 2/177 Orlando St, Coffs Harbour NSW 2450

Tel: (02) 6651 5333

BALLINA
18 Endeavour CI,
Ballina NSW 2478
Tel: (02) 6686 6644
Fax: (02) 6686 9383

QUEENSLAND

BRISBANE SOUTH
30 Dulacca St,
Acacia Ridge QLD, 4110
Tel: (07) 3086 2900
Fax: (07) 3345 7961

BRISBANE NORTH 354 Bilsen Rd, Geebung OLD 4034 Tel: (07) 3086 2950 Fax: (07) 3865 7896

SUNSHINE COAST
Unit 2, 35 Premier Circuit
Warana, OLD, 4575
Tel: (07) 3086 2980
Fax: (07) 3345 7961

HERVEY BAY Shop 3, 109 Beach Road Pialba, QLD, 4655 Tel: (07) 4194 5019

GOLD COAST
Unit 1/ 48 Olympic Circuit,
Southport QLD 4215
Tel: (07) 5532 6070
Fax: (07) 5503 1803

TOOWOOMBA

1/259 James St,
Toowoomba QLD 4350
Tel: (07) 3086 2990
Fax: (07) 4632 6473

SOUTH AUSTRALIA

ADELAIDE 21-23 Tikalara St, Regency Park SA 5010 Tel: (08) 8118 5600 Fax: (08) 8346 3236

TASMANIA

18 Mertonvale Circuit, Kingston TAS 7050 Tel: (03) 6229 4733 Fax: (03) 6229 3676

Unit 8/34 Innocent St, King Meadows TAS 7250 Tel: (03) 6334 8922 Fax: (03) 6229 3676

National Distribution Network



CALL 1300 133 120

Call this number to be automatically diverted to your local Aidacare office (for the price of a local call).

DVA Podiatry Product List



Product Index

DVA: AB06 **Bed Cradles** P 4 DVA: AB11 Sheepskin Rugs / Foot / Heel / Elbow Pads DVA: AB14 P 8 **Pressure Care Mattresses** DVA: AB17 **Heel Elevators for Pressure Care** P 13 DVA: AJ06 Footwear (Including Cast Boots / Shoes) P 14 DVA: AR04 Foot Orthoses / Orthotics

DVA: AR19 **Ankle Foot Orthoses (AFO)** DVA: AU10 Long Handled Toe Wiper

DVA: AU13 P 19 **Donning / Doffing Aids**

DVA: AU14







• Powder Coated Mild Steel Construction



aspire



| Code | Description | Overall Height | Overall Width | Overall Depth | DVA: |
|-----------|--------------------------|----------------|---------------|---------------|------|
| BEA014605 | Bed Cradle White (Fixed) | 490 mm | 470 mm | 320 mm | AB06 |



Bed Cradle - Adjustable

- Height Adjustable
- Powder Coated Mild Steel Construction

| Code | Description | Overall Height | Overall Width | Overall Depth | DVA: |
|-----------|-------------------------------|----------------|---------------|---------------|------|
| BEA014610 | Bed Cradle White (Adjustable) | 440mm - 640mm | 370 mm | 700 mm | AB06 |

eRAP: DVA Prescriber Portal

eRAP is Aidacare's mobile platform for Health Providers to submit RAP direct order forms. Prescribers can access our eRAP via the web portal or via our iPad app.



For more information, visit erap.aidacare.com.au or call 1300 888 052

Sheepskin Rugs / Foot / Heel / Elbow Pads







Sheepskin Leg Rest Cover

- High quality medical grade sheepskin covers provide natural comfort for any wheelchair
- Provide excellent skin protection and comfort

| Code | Description | DVA: |
|-----------|----------------------------|------|
| PTA522100 | Wheelchair Footplate Cover | AB11 |
| PTA522101 | Wheelchair Calf Protector | AB11 |



Sheepskin Hi Temp - Cushion Seat Back

- The Cushion It Sheepskin Cover can be used directly on a seat or back canvas of a wheelchair
- Can be placed over standard foam cushions to increase pressure care and absorb perspiration

| Code | Description | Length | Width | DVA: |
|-----------|--------------------------------------|--------|--------|------|
| PTA521300 | Sheepskin HiTemp - Cushion Seat Back | 500 mm | 500 mm | AB11 |



Sheepskin Hi Temp Medical Overlay

- Medical grade sheepskin is a versatile product that can be used in the bed, on a chair or on a car seat
- · Assists in pressure reduction, will absorb water, transfer heat away from the body and reduce shear and friction

| Code | Description | Length | Width | DVA: |
|-----------|---------------------------|---------|--------|------|
| PTA501500 | Medical Sheepskin Overlay | 1000 mm | 650 mm | AB11 |



Sheepskin Natural Gold Overlay

- Genuine Australian lambswool overlay provides additional comfort to an existing mattress
- Soft pile allows air circulation to keep skin dry

| Code | Description | Length | Width | DVA: |
|-----------|-------------------|---------|--------|------|
| PTA502500 | Sheepskin Overlay | 1000 mm | 750 mm | AB11 |



Sheepskin Hi Temp Classic Mattress Overlay

- A full covered medical sheepskin bed overlay, designed to cover head to toe
- Suitable for long term patients where shear and friction are a concern

| Code | Description | Length | Width | DVA: |
|-----------|----------------------------|---------|--------|------|
| PTA521400 | Sheepskin Mattress Overlay | 1250 mm | 750 mm | AB11 |

Sheepskin Rugs / Foot / Heel / Elbow Pads







Pressure Care Assist

- Can be used anywhere in the bed to assist with pressure care
- Made from medical grade sheepskin
- Helps reduce the risk of pressure area breakouts in the heels, sacrum, shoulders or head

| Code | Description | DVA: |
|-----------|----------------------|------|
| PTA521600 | Pressure Care Assist | AB11 |



Aspire Heel Protector



- Intelligent leg support zone offloads pressure from the heel and redistributes it to the larger surface area of the lower leg.
- Soft breathable natural cotton fabric outer.
- Cushioning silicone fibre inner promotes airflow for skin microclimate management.

| Code | Description | Quantity | Sizing | DVA: |
|-----------|-----------------------|----------|---|------|
| PTA503360 | Aspire Heel Protector | 1 x Pair | Universal design fits majority of users | AB11 |



Sheepskin Short Slipper Boot

- Assists in the prevention of pressure ulcers on the feet with an existing injury or condition
- Designed for use primarily in bed, wheelchair use or recliner chair use

| Code | Description | Shoe Size US | DVA: |
|-----------|-------------------------------|--------------|------|
| PTA521897 | Short Slipper Boot - XX-Small | 3 - 4 | AB11 |
| PTA521898 | Short Slipper Boot - X-Small | 5 - 6 | AB11 |
| PTA521899 | Short Slipper Boot - Small | 7 - 8 | AB11 |
| PTA521900 | Short Slipper Boot - Medium | 9 - 10 | AB11 |
| PTA521901 | Short Slipper Boot - Large | 11 - 12 | AB11 |
| PTA521902 | Short Slipper Boot - X-Large | 13 - 14 | AB11 |
| PTA521903 | Short Slipper Boot - XX-Large | 14 - 15 | AB11 |



Sheepskin Wrap Around Boot

- Wrap Around Sheepskin Boots designed for those at risk of foot ulcers and other types of tissue damage
- Constructed with low profile EVA sole, these boots can be worn outside

| Code | Description | Shoe Size US | DVA: |
|-----------|-----------------------------|--------------|------|
| PTA521704 | Wrap Around Boot - XX-Small | 3 - 4 | AB11 |
| PTA521705 | Wrap Around Boot - X-Small | 5 - 6 | AB11 |
| PTA521706 | Wrap Around Boot - Small | 7 - 8 | AB11 |
| PTA521707 | Wrap Around Boot - Medium | 9 - 10 | AB11 |
| PTA521708 | Wrap Around Boot - Large | 11 - 12 | AB11 |
| PTA521709 | Wrap Around Boot - X-Large | 13 - 14 | AB11 |
| PTA521710 | Wrap Around Boot - XX-Large | 14 - 15 | AB11 |



Sheepskin Pressure Care Boot

- Sheepskin Boots designed for those at risk of foot ulcers and other types of tissue damage
- Constructed with low profile EVA sole, these boots can be worn outside

| Code | Description | Shoe Size US | DVA: |
|-----------|-------------------------------|--------------|------|
| PTA521697 | Pressure Care Boot - XX-Small | 3 - 4 | AB11 |
| PTA521698 | Pressure Care Boot - X-Small | 5 - 6 | AB11 |
| PTA521699 | Pressure Care Boot - Small | 7 - 8 | AB11 |
| PTA521700 | Pressure Care Boot - Medium | 9 - 10 | AB11 |
| PTA521701 | Pressure Care Boot - Large | 11 - 12 | AB11 |
| PTA521702 | Pressure Care Boot - X-Large | 13 - 14 | AB11 |
| PTA521703 | Pressure Care Boot - XX-Large | 14 - 15 | AB11 |

Sheepskin Rugs / Foot / Heel / Elbow Pads







Sheepskin Tall Slipper Boot

- This product has been designed to protect the entire lower limb
- The Tall Slipper Sheepskin Boot can be worn in bed, in a chair or wheelchair and allows the use of compression bandages or support hosiery

| Code | Description | DVA: |
|-----------|---------------------------------------|------|
| PTA520500 | Tall Slipper Boot - One Size Fits All | AB11 |



Sheepskin Sovereign Snug

- Beautifully crafted Sheepskin Snugs will accommodate mild foot oedema whilst offering support and comfort
- Made especially from carefully selected lambskin material

| Code | Description | Shoe Size US | DVA: |
|-----------|-------------------------|--------------|------|
| PTA521905 | Sovereign Snug - Small | 7 - 8 | AB11 |
| PTA521906 | Sovereign Snug - Medium | 9 - 10 | AB11 |
| PTA521907 | Sovereign Snug - Large | 11 - 12 | AB11 |



Sheepskin Heel Supports

- Protect the heel from ulcerations caused by bumps and knocks
- Made from medical grade sheepskin
- Wraps around the heel and ankle leaving the rest of the foot unencumbered

| Code | Description | Show Size US | DVA: |
|-----------|-------------------------------|--------------|------|
| PTA521199 | Heel Supports - Small (Pair) | 7 - 8 | AB11 |
| PTA521200 | Heel Supports - Medium (Pair) | 9 - 10 | AB11 |
| PTA521201 | Heel Supports - Large (Pair) | 11 - 12 | AB11 |



Sheepskin Elbow Protector

- Provides elbow pressure point control
- Assists in avoiding skin tears, bruising and other tissue damage
- Easy and convenient velcro fastening
- Sold as a pair

| Code | Description | DVA: |
|-----------|-------------------------------------|------|
| PTA521100 | Elbow Protector (One Size Fits All) | AB11 |



Sheepskin Palm Protector

- Provides comfort for patients with hand contractures
- Designed to offer moisture absorption, reducing the likelihood of maceration
- One size fits all

| Code | Description | DVA: |
|-----------|-------------------------------|------|
| PTA521104 | Palm Protector (Ambidextrous) | AB11 |









Lifecomfort Pressure Safety Underlay

- Provides "bottoming out" protection for air alternating mattress replacements
- Minimal height addition to overall pressure care system
- Outer cover is antimicrobial, anti bacterial, water resistant and fire retardant
- Four way rotation ensures longevity of foam performance
- Vapour permeable (breathable)

| Code | Description | Width | Length | Height | DVA: |
|-----------|--|---------|---------|--------|------|
| BEM047050 | Lifecomfort Pressure Safety Underlay - Single | 900 mm | 2000 mm | 50 mm | AB14 |
| BEM047060 | Lifecomfort Pressure Safety Underlay - King Single | 1050 mm | 2030 mm | 50 mm | AB14 |



M/ Overlay Naturelle Latex

- · Anti-microbial and dust-mite resistant to alleviate allergies
- Hypoallergenic, hygienic, and antibacterial
- Optimal orthopaedic support (naturally conforming to body contours)
- Superior pressure relief that helps bring oxygen and nutrients to muscles

| Code | Description | Width | Length | Height | S.W.L | DVA: |
|-----------|-------------------------------------|---------|---------|--------|--------|------|
| BEM047181 | M/ Overlay Naturelle Latex - Single | 905 mm | 1880 mm | 50 mm | 100 kg | AB14 |
| BEM047182 | M/ Overlay Naturelle Latex - Double | 1355 mm | 1880 mm | 50 mm | 100 kg | AB14 |
| BEM047183 | M/ Overlay Naturelle Latex - Queen | 1505 mm | 2000 mm | 50 mm | 100 kg | AB14 |
| BEM047184 | M/ Overlay Naturelle Latex - King | 1805 mm | 2000 mm | 50 mm | 100 kg | AB14 |



Mattress Overlay Spenco

- Silcone coated fibre filled overlay
- Cylindrical hollow core fibres provide greater comfort and shear reduction
- Quality breathable cotton cover
- Heavy duty stitching
- · Supplied with heat sealed waterproof cover

| Code | Description | Width | Length | Warranty | S.W.L | DVA: |
|-----------|-------------------------|--------|---------|----------|--------|------|
| PTM636160 | Mattress Overlay Spenco | 910 mm | 1980 mm | 1 year | 100 kg | AB14 |



Mattress Overlay Silicone Fibre

- Fits under bed sheet to reduce friction and shear for the whole body
- Silicone coated cylindrical hollow core fibres construction
- · Assists with geriatric and palliative care as silicone fibres do not absorb perspiration or urine
- · Helps improve sleep by adding extra softness and support to the bed, even if lying in one position for extended periods

| Code | Description | Width | Length | Height | S.W.L | DVA: |
|-----------|---------------------------------|--------|---------|--------|--------|------|
| PTM636170 | Mattress Overlay Silicone Fibre | 910 mm | 1900 mm | 100 mm | 100 kg | AB14 |



Convoluted Overlay

- Fits under bed sheet to re-distribute pressure points for therapeutic pressure relief and comfort
- Improves air flow around the mattress helping to maintain a constant body temperature and improve sleep quality
- Egg shell shape/convoluted foam

| Code | Description | Width | Length | Height | Warranty | DVA: |
|-----------|-----------------------------|---------|---------|--------|----------|------|
| PTM636260 | Convoluted Overlay - Single | 915 mm | 1880 mm | 50 mm | 1 year | AB14 |
| PTM636270 | Convoluted Overlay - Double | 1370 mm | 1880 mm | 50 mm | 1 year | AB14 |
| PTM636280 | Convoluted Overlay - Queen | 1520 mm | 2030 mm | 50 mm | 1 year | AB14 |

Pressure Care Mattresses









Lifecomfort Foam Mattress

- Pressure care mattress
- Two layer foam support design for comfort and support
- Full size memory foam overlay
- Two way waterproof stretch cover
- Vinyl bottom cover

| Code | Description | Width | Length | Height | SWL | DVA: |
|-----------|--------------------------------|---------|---------|--------|--------|------|
| BEM046410 | Lifecomfort - AC Series Single | 900 mm | 2000 mm | 130 mm | 120 kg | AB14 |
| BEM046425 | Lifecomfort - Small | 820 mm | 1920 mm | 130 mm | 120 kg | AB14 |
| BEM046460 | Lifecomfort - King Single | 1050 mm | 2030 mm | 130 mm | 120 kg | AB14 |



Lifecomfort Deluxe Foam Mattress



- Three layer foam design for enhanced comfort and support
- Full size memory foam overlay
- Two way waterproof stretch cover
- Vinyl bottom cover

| Code | Description | Width | Length | Height | SWL | DVA: |
|-----------|---------------------------------------|---------|---------|--------|--------|------|
| BEM046510 | Lifecomfort Deluxe - AC Series Single | 900 mm | 2000 mm | 150 mm | 150 kg | AB14 |
| BEM046525 | Lifecomfort Deluxe - Small | 820 mm | 1920 mm | 150 mm | 150 kg | AB14 |
| BEM046455 | Lifecomfort Deluxe - King Single | 1050 mm | 2030 mm | 150 mm | 150 kg | AB14 |




Lifecomfort Bariatric Foam Mattress



- Pressure reducing bariatric mattress
- Extra weight capacity up to 300 kg
- Four layer foam support 200 mm high for bariatric use with maximum comfort and support
- Full size memory foam overlay
- Foam side walls for extra stability on the edge of the mattress

| Code | Description | Width | Length | Height | SWL | DVA: |
|-----------|-----------------------|---------|---------|--------|--------|------|
| BEM046700 | Lifecomfort Bariatric | 1050 mm | 2030 mm | 200 mm | 300 kg | AB14 |



Lifecomfort Pressure Reducing Foam Mattress



- Pressure reducing mattress
- Four layer foam support design for maximum comfort and support
- Full size memory foam overlay
- Foam side walls for extra stability on the edge of the mattress
- Two way waterproof stretch cover with welded seams

| Code | Description | Width | Length | Height | SWL | DVA: |
|-----------|---|---------|---------|--------|--------|------|
| BEM046610 | Lifecomfort Pressure Reducing - AC Series Single | 900 mm | 2000 mm | 150 mm | 150 kg | AB14 |
| BEM046660 | Lifecomfort Pressure Reducing - King Single | 1050 mm | 2030 mm | 150 mm | 150 kg | AB14 |




Lifecomfort Foam Mattress with Safety Sides



- Pressure care mattress with sidewall support to help minimise the risk of falls from bed
- Two layer foam support design for comfort and support
- Full size memory foam overlay
- Two way waterproof stretch cover
- · Vinyl bottom cover

| Code | Description | Width | Length | Height | Safety Sides | SWL | DVA: |
|-----------|-------------------------------|--------|---------|--------|--------------|--------|------|
| BEM046475 | Lifecomfort with Safety Sides | 900 mm | 1940 mm | 130 mm | 220 mm | 120 kg | AB14 |









Lifecomfort Foam Mattress Extension Bolster

🧊 aspire

· Helps your facility deal with the variability of bed lengths reducing the risk of injury to the patient.

| Code | Description | Width | Length | Height | DVA: |
|-----------|--|---------|--------|--------|------|
| BEM046762 | Mattress Extension Bolster - Short - Single | 900 mm | 100 mm | 130 mm | AB14 |
| BEM046764 | Mattress Extension Bolster - Short - King Single | 1050 mm | 100 mm | 130 mm | AB14 |
| BEM046772 | Mattress Extension Bolster - Long - Single | 900 mm | 200 mm | 130 mm | AB14 |
| BEM046774 | Mattress Extension Bolster - Long - King Single | 1050 mm | 200 mm | 130 mm | AB14 |



Lifecomfort Hybrid Mattress Replacement

👣 aspire

- Intelligent Hose Management
- Memory Foam Top Layer
- Optimised Heel Protection
- Integrated Shear Reduction

| Code | Description | Width | Length | Height | No. Of Cells | SWL |
|-----------|---|---------|---------|--------|--------------|--------|
| PTM638010 | Hybrid Mattress Replacement - Single | 895 mm | 1980 mm | 170 mm | 9 | 180 kg |
| PTM638015 | Hybrid Mattress Replacement - King Single | 1050 mm | 1980 mm | 170 mm | 9 | 180 kg |



Lifecomfort Active Air 4 Overlay



- The Active Air 4 is a lightweight, easy to use Overlay that is placed on top of an existing foam mattress. Featuring COMFORT CONTROL and offering a suite of other features, the Active Air 4 caters to users considered to be "At Risk" of Pressure Injury.
- An ideal solution for Home Care or Facility environments.

| Code | Description | Width | Length | Height | Weight | No. of Cells | S.W.L |
|-----------|----------------------|--------|---------|--------|--------|--------------|--------|
| PTM639100 | Active Air 4 Overlay | 850 mm | 2000 mm | 100 mm | 4.9 kg | 19 | 135 kg |



Lifecomfort Active Air 5 Overlay



- The Active Air 5 provides an effective solution for higher weight users considered to be "At Risk" of Pressure Injury.
- Featuring AUTO DETECT and INTELLIGENT STATIC MODE which provide lower interface pressure and enhanced comfort, the Active Air 5 is suitable for use in a wide range of care settings from Home Care and Acute to Long Term care.

| Code | Description | Width | Length | Height | Weight | No. of Cells | S.W.L |
|-----------|------------------------------------|---------|---------|--------|--------|--------------|--------|
| PTM639200 | Active Air 5 Overlay - Single | 850 mm | 2000 mm | 130 mm | 5.3 kg | 16 | 175 kg |
| PTM639210 | Active Air 5 Overlay - King Single | 1050 mm | 2000 mm | 130 mm | 5.7 kg | 16 | 175 kg |



Lifecomfort Active Air 8 Mattress Replacement



- The Active Air 8 is a dynamic mattress replacement suitable for users at "High Risk" of Pressure Injury.
- Featuring AUTODETECT, and AUTOPROTECT, together with CELL-IN-CELL technology, the Active Air 8 provides the highest levels of care into Acute and Long Term care environments. The Active Air 8's innovative features ensure lowest interface pressures with an easy-to-use digital pump.

| Code | Description | Width | Length | Height | Weight | No. of Cells | S.W.L |
|-----------|--|---------|---------|--------|---------|--------------|--------|
| PTM639300 | Active Air 8 Mattress Replacement - Single | 900 mm | 2000 mm | 200 mm | 9.5 kg | 20 | 200 kg |
| PTM639310 | Active Air 8 Mattress Replacement - King Single | 1050 mm | 2000 mm | 200 mm | 10.9 kg | 20 | 200 kg |

Pressure Care Mattresses







ROHO DRY FLOATATION® Mattress Overlay

- Dry flotation mattress system allows user to lie or sit comfortably for long periods
- Interconnected and flexible air cells provide protection and relief for user at risk of skin breakdown or already with early stage ulcers
- Adjustable air pressure to suit the user
- · Flexibility provides effective friction and shear relief

| Code | Description | Height | Width | Length | SWL | Warranty | DVA: |
|-----------|---|--------|--------|---------|-----------|----------|------|
| PTM636335 | ROHO DRY FLOATATION® Mattress Overlay (4 sections) | 80 mm | 860 mm | 1960 mm | Unlimited | 2 Years | AB14 |
| PTM636320 | ROHO DRY FLOATATION® Mattress Section (1 section) | 80 mm | 860 mm | 490 mm | Unlimited | 2 Years | AB14 |



ROHO® Levelling Pad

- Dry flotation mattress system allows user to lie or sit comfortably for long periods
- Interconnected and flexible air cells provide protection and relief for user at risk of skin breakdown or already with early stage ulcers
- Adjustable air pressure to suit the user
- · Flexibility provides effective friction and shear relief

| Code | Description | Width | Length | Height | SWL | Warranty | DVA: |
|-----------|-----------------------------|--------|---------|--------|--------|----------|------|
| PTM636360 | ROHO® Levelling Pad - Long | 864 mm | 2032 mm | 80 mm | 100 kg | 2 Years | AB14 |
| PTM636220 | ROHO® Levelling Pad - Short | 863 mm | 508 mm | 80 mm | 100 kg | 2 Years | AB14 |



Air-Lo 3" (75mm) Overlay

- Static head section eliminates patient dizziness and disorientation
- 21 pure TPU transverse air cells alternate to provide maximum comfort prevention and therapy
- High quality two way stretch top cover is fully welded and can be high temperature washed for total infection control

| Code | Description | Width | Length | Height | CycleTime | SWL | Warranty | DVA: |
|-----------|---|--------|---------|--------|-----------|--------|----------|------|
| PTM636629 | Air-Lo 3 | 900 mm | 2000 mm | 75 mm | 14 mins | 120 kg | 2yr/1yr | AB14 |
| PTM636625 | AirLo -Thermal Hyperstretch Microfibre Coverlet | | | | | | | |



Air-Lo 5" (125mm) Overlay

- Static head section eliminates patient dizziness and disorientation & releases with CPR
- 18 pure TPU transverse air cells alternate to provide maximum comfort prevention and therapy
- High quality two way stretch top cover is fully welded and can be high temperature washed for total infection control
- All waterfalls are turned and welded for maximum infection control and zip protection

| Code | Description | Width | Length | Height | CycleTime | SWL | Warranty | DVA: |
|-----------|---|--------|---------|--------|-----------|--------|----------|------|
| PTM636630 | Air-Lo 5 | 900 mm | 2000 mm | 125 mm | 14 mins | 150 kg | 2yr/1yr | AB14 |
| PTM636625 | AirLo -Thermal Hyperstretch Microfibre Coverlet | | | | | | | AB14 |



Air-Lo 8" (200mm) Mattress Replacement

- Static head section eliminates patient dizziness and disorientation
- 20 pure TPU transverse air cells alternate to provide maximum comfort prevention and therapy
- High quality two way stretch top cover is fully welded and can be high temperature washed for total infection control
- All waterfalls are turned and welded for maximum infection control and zip protection

| Code | Description | Width | Length | Height | CycleTime | SWL | Warranty | DVA: |
|---|-------------|--------|---------|--------|-----------|--------|----------|------|
| PTM636635 | Air-Lo 8 | 900 mm | 2000 mm | 200 mm | 14 mins | 200 kg | 2yr/1yr | AB14 |
| PTM636625 AirLo -Thermal Hyperstretch Microfibre Coverlet | | | | | | AB14 | | |







Ultimate 75 (3") Overlay

- A thinner mattress overlay support system for patients at medium risk of pressure ulcers
- An economic and comfortable solution to pressure injury prevention
- Lightweight, manoeuvrable, versatile and easy to place over any existing mattress
- Supports users to 120 kg
- End straps securely fasten the system to the mattress

| Code | Description | Width | Length | Height | CycleTime | SWL | Warranty | DVA: |
|-----------|-------------|--------|---------|--------|-------------|--------|----------|------|
| PTM636560 | Ultimate 75 | 880 mm | 2000 mm | 75 mm | 8 - 10 mins | 120 kg | 2yr/1yr | AB14 |



Ultimate 125 (5") Overlay

- · A medium mattress overlay support system for patients at medium risk of developing pressure ulcers
- 2 cell cycle 21 18 cells
- Prevention of pressure injury is the core strategy of this product
- Digital control unit with audio alarm to alert carer to malfunction
- · Static button to activate full inflation of all cells enabling ease of process when transferring and treating patients

| Code | Description | Width | Length | Height | CycleTime | SWL | Warranty | DVA: |
|-----------|----------------------|--------|---------|--------|-----------|--------|----------|------|
| PTM636530 | Ultimate 125 Digital | 880 mm | 2000 mm | 125 mm | 12 mins | 150 kg | 2yr/1yr | AB14 |



Ultimate 200 (8") Mattress Replacement

- An economic, full mattress support system for patients at high risk of developing pressure injuries
- 2 cell cycle 21 cells
- Increased cell depth is ideal for the heavier patient, prevents "bottoming out" and enhances therapeutic
- · Heat sealed seams, antibacterial taping and waterfall skirts facilitate effective infection control and prevention

| Code | Description | Width | Length | Height | CycleTime | SWL | Warranty | DVA: |
|-----------|----------------------|--------|---------|--------|-----------|--------|----------|------|
| PTM636540 | Ultimate 200 Digital | 880 mm | 2000 mm | 200 mm | 8-10 mins | 180 kg | 2yr/1yr | AB14 |



Ultimate 225 (9") Mattress Replacement

- A professional, clinical care, full mattress support system for patients at very high risk of developing pressure injuries
- 3 cell cycle 21 cells
- Unique heel zone ensures pressure alternation of high risk areas of the lower leg
- · Lateral side bolsters for additional patient stability
- Excellent healing surface for patients receiving treatment for existing wounds

| Code | Description | Width | Length | Height | CycleTime | SWL | Warranty | DVA: |
|-----------|-------------------------|--------|---------|--------|------------|--------|----------|------|
| PTM636545 | Ultimate 225 Digital | 880 mm | 2000 mm | 225 mm | 8 -10 mins | 180 kg | 2yr/1yr | AB14 |

Heel Elevators for Pressure Care







Leg Relaxer

- Inserted into bed to assist users body positioning and posture
- · Aids circulation in the lower limbs and provides even weight distribution alleviating lower back strain
- Ideal for those with varicose veins or swollen legs, and during pregnancy

| Code | Description | Width | Length | Height | DVA: |
|-----------|----------------------------------|--------|--------|--------|------|
| BEP056200 | Leg Relaxer with Quilted Slip | 450 mm | 710 mm | 200 mm | AB17 |
| BEP056250 | Leg Relaxer with Waterproof Slip | 450 mm | 710 mm | 200 mm | AB17 |



Bed Heel Elevator

- Wedge-shaped contoured foam provides excellent support for the lower leg
- Rigid foam aids blood flow and reduces pressure

| Code | Description | DVA: |
|-----------|-------------------|------|
| PTA509050 | Bed Heel Elevator | AB17 |



ROHO Heelpad

- ROHO Heelpad offers an adjustable therapeutic environment which aids in the healing of heel and elbow ischemic ulcers
- Made from soft, latex free neoprene the heelpad is extremely durable

| Code | Description | DVA: |
|-----------|--------------|------|
| PTC541600 | ROHO Heelpad | AB17 |

eRAP: DVA Prescriber Portal

eRAP is Aidacare's mobile platform for Health Providers to submit RAP direct order forms. Prescribers can access our eRAP via the web portal or via our iPad app.



For more information, visit erap.aidacare.com.au or call 1300 888 052



Footwear (Including Cast Boots / Shoes)







Darco All Purpose Boot

- Rocker Sole reduces plantar pressure under forefoot and heel by over 25%
- Square Toe Design acts as a bumper to protect the toe
- Strap seats foot firmly in shoe to reduce heel slippage and eliminate friction
- Universal Left/Right fit

| Codo | Ciro | Shoe | DV/A. | |
|-----------|------|-----------|-----------|------|
| Code | Size | Female | Male | DVA: |
| FBS280240 | XS | 4 - 7 | - | AJ06 |
| FBS280241 | S | 7.5 - 10 | 6 - 8 | AJ06 |
| FBS280242 | М | 10.5 - 13 | 8.5 - 10 | AJ06 |
| FBS280245 | L | 13+ | 10.5 - 12 | AJ06 |
| FBS280246 | XL | - | 12.5 - 14 | AJ06 |

Note: Sizing by US shoe size.



Darco Heelwedge Shoe

- Effective off-loading Footwear
- Reduces Heel pressure by over 26%
- Square Toe Design acts as a bumper to protect the toes and provides better fit
- Ankle Strap seats foot firmly in shoe to reduce heel slippage and eliminate friction
- Zoned Outsole has an aggressive tread under the mid-foot where traction is needed most

| Cada | Ciro | Shoe | DV/A. | |
|-------------|------|-----------|-----------|------|
| Code | Size | Female | Male | DVA: |
| FBS301577 | XS | 4 - 7 | - | AJ06 |
| FBS3015771B | S | 7.5 - 10 | 6 - 8 | AJ06 |
| FBS3015772B | M | 10.5 - 13 | 8.5 - 10 | AJ06 |
| FBS3015773B | L | 13+ | 10.5 - 12 | AJ06 |
| FBS3015774B | XL | - | 12.5 - 14 | AJ06 |

Note: Sizing by US shoe size.



Darco Med Surg Shoe

- Square Toe Design acts as a bumper to provide extra protection and facilitates better left/right fit
- MetaShank™ protection provides more rigid control under the metatarsal heads
- Rocker Sole reduces plantar pressure on the heel and forefoot by over 25%
- Fully compatible with DARCO's Peg-Assist™ off-loading insole

| | Code | | Size | Shoe Size | | |
|-------------|-------------|------|------|-----------|-----------|--|
| Female | Female Male | | Size | Female | Male | |
| FBS301579 | FBS301578 | AJ06 | S | 4 - 6 | 6 - 8 | |
| FBS3015792B | FBS3015782B | AJ06 | M | 6.5 - 8 | 8.5 - 10 | |
| FBS3015793B | FBS3015783B | AJ06 | L | 8.5 - 10 | 10.5 - 12 | |
| - | FBS3015784B | AJ06 | XL | - | 12.5 - 14 | |

Note: Sizing by US shoe size.



Darco Orthowedge Shoe

- 15° Wedge Sole effectively shifts body weight to the midfoot and heel
- Reduces forefoot pressure by as much as 57%
- Square Toe Design acts as a bumper to protect the toes when k wires are present and provides better left/ right fit

| 0.1 | C: | Shoe Size | | | |
|-------------|------|-----------|-----------|------|--|
| Code | Size | Female | Male | DVA: | |
| FBS301580 | XS | 4-7 | - | AJ06 | |
| FBS3015801B | S | 7.5-10 | 6 - 8 | AJ06 | |
| FBS3015802B | M | 10.5-13 | 8.5 - 10 | AJ06 | |
| FBS3015803B | L | 13÷ | 10.5 - 12 | AJ06 | |
| FBS3015804B | XL | - | 12.5 - 14 | AJ06 | |

Note: Sizing by US shoe size.

Footwear (Including Cast Boots / Shoes)







Darco Soffie Shoe

- Economical, post-op shoe indicated for soft tissue trauma or surgery
- Semi-rigid metatarsal shank provides control of forefoot motion
- Tissue friendly upper material conforms to bony abnormalities and is breathable
- Easy to fasten hook and loop straps
- Female/Male specific for optimal fit

| Co | Code | | Shoe Size | | |
|-----------|-----------|------|-----------|-----------|------|
| Female | Male | Size | Female | Male | DVA: |
| FBS280230 | FBS280233 | S | 4 - 6 | 6 - 8 | AJ06 |
| FBS280231 | FBS280234 | M | 6.5 - 8 | 8.5 - 10 | AJ06 |
| FBS280232 | FBS280235 | L | 8.5 - 10 | 10.5 - 12 | AJ06 |
| - | FBS280236 | XL | - | 12.5 - 14 | AJ06 |

Note: Sizing by US shoe size.



Cosmac Cast Shoe

- Curved 'rocker' sole for easy ambulation
- · Canvas upper with velcro closures at front and back
- Symmetrical shape fits left & right foot

| Code | Description | Size | Shoe Length (cm) | DVA: |
|------------|------------------|------|---------------------|------|
| FBS301583S | Cosmac Cast Shoe | S | 23 | AJ06 |
| FBS301583M | Cosmac Cast Shoe | M | 27 | AJ06 |
| FBS301583L | Cosmac Cast Shoe | L | 32 | AJ06 |

Note: Sizing by US shoe size.



Darco Body Armor Cast Shoe

- Ultimate cast protection in a new modern design that will appeal to all patients
- · Soft but durable EVA provides shock absorption while protecting the cast from the elements
- · Unique bungee closure system is easy to use, easy to adjust and prevents the heel from slipping out of the cast shoe
- Soft, durable TPR

| Code Size | | Shoe | DVA: | |
|------------|------|-----------|-----------|------|
| Code | Size | Female | Male | DVA: |
| FBS301582 | XS | Up to 5 | - | AJ06 |
| FBS3015821 | S | 5.5 - 8 | - | AJ06 |
| FBS3015822 | M | 8.5 - 10 | 6 - 8 | AJ06 |
| FBS3015823 | L | 10.5 - 13 | 8.5 - 11 | AJ06 |
| FBS3015824 | XL | 13.5+ | 11.5 - 14 | AJ06 |

Note: Sizing by US shoe size.



Darco Slimline Cast Boot

- Be adapted to fiberglass casts, bulky compression bandages and a variety of woundcare applications, making it one of the most versatile cast shoes on the market
- Forefoot Closure can be adjusted to accommodate casts or bandages
- Symmetrical shape fits left & right foot
- Higher Ankle Strap adds stability
- Rocker Sole allows for a smooth roll throughout the gait pattern

| Code | Size | Shoe Size | | | |
|-------------|------|-----------|-----------|------|--|
| Code | Size | Female | Male | DVA: | |
| FBS301583 | XS | Up to 5 | - | AJ06 | |
| FBS3015831B | S | 5.5 - 8 | - | AJ06 | |
| FBS3015832B | M | 8.5 - 10 | 6 - 8 | AJ06 | |
| FBS3015833B | L | 10.5 - 13 | 8.5 - 11 | AJ06 | |
| FBS3015834B | XL | 13.5+ | 11.5 - 14 | AJ06 | |

Note: Sizing by US shoe size.







BodyAssist Heel Spur Heel Cups

- Silicone Medical Gel provide exceptional heel support, protection & comfort
- Ideal for sports use, gardening, factory staff, construction workers, office staff & for people who stand on hard surfaces
- Heel Spur Pads are made from a unique Medical Silicone Gel compound that assists to dissipate shock & shear, caused by impact to the heel & foot
- They emulate the fatty tissue, usually found on the underfoot
- Damaging stress forces normally absorbed by the body's musculo-skeletal system, are now absorbed by the pad & away from the body

| Code | Size | Shoe Size | DVA: |
|------------|------|-----------|------|
| FBS281210 | S | 37- 40 | AR04 |
| FBS2812102 | M | 40.5 - 42 | AR04 |
| FBS2812103 | L | - | AR04 |



Valgus 4050 Arch/Metatarsal Support

- Anatomically designed moulded rubber with leather cover
- Metatarsal dome
- Available in narrow and wide fittings
- Supplied in pairs

| Shoe Size | | Code | |
|-------------|------------------|------------------|------|
| Si lue Size | Narrow | Wide | DVA: |
| 33 | FBS301590 | - | AR04 |
| 34 | FBS301590/100/34 | - | AR04 |
| 35 | FBS301590/100/35 | - | AR04 |
| 36 | FBS301590/100/36 | - | AR04 |
| 37 | FBS301590/100/37 | - | AR04 |
| 38 | FBS301590/100/38 | FBS301590/200/38 | AR04 |
| 39 | FBS301590/100/39 | FBS301590/200/39 | AR04 |
| 40 | FBS301590/100/40 | FBS301590/200/40 | AR04 |
| 41 | FBS301590/100/41 | FBS301590/200/41 | AR04 |
| 42 | FBS301590/100/42 | FBS301590/200/42 | AR04 |
| 43 | - | FBS301590/200/43 | AR04 |
| 44 | - | FBS301590/200/44 | AR04 |
| 45 | - | FBS301590/200/45 | AR04 |
| 46 | - | FBS301590/200/46 | AR04 |



BodyAssist Heel Spur Pads

- These Heel Spur Pads gently cushion the heel, while dissipating shock and shear generated by heel strike
- A removable 'tear drop' central plug allows the heel centre to be suspended around a supportive ring of silicone gel
- The pads fit snugly into the heel of a shoe
- Shock attenuation properties of this material have been proven effective in modifying heel strike trauma

| Code | Size | Shoe Size | DVA: |
|------------|------|-----------|------|
| FBS301598 | S | Up to 38 | AR04 |
| FBS3015982 | M | 38.5 - 42 | AR04 |
| FBS3015983 | L | 42.5+ | AR04 |

Foot Orthoses / Orthotics







BodyAssist Gel Shoe Insole

- · Provides overall support, protection, comfort & help reduce painful shock waves to the body & lighten the load on joints
- Ideal for sports use, gardening, factory staff, construction workers, office staff & for people who stand on hard surfaces
- · Made from a unique Medical Silicone Gel compound that assists to dissipate shock & shear, caused by impact to the heel & foot
- Suitable for managing minor cases of splayfoot

| Code | Size | Shoe Size | DVA: |
|-----------|------|-----------|------|
| FBS281200 | XS | 37 - 38 | AR04 |
| FBS281201 | S | 39 - 41 | AR04 |
| FBS281202 | M | 41 - 43 | AR04 |
| FBS281203 | L | 44 - 46 | AR04 |
| FBS281204 | XL | 47+ | AR04 |



Trio Heel Elevators

- Trio elevators come in different heights and should be placed under the heel area of the orthotic
- Choose between sizes small, medium or large
- They are manufactured from extra firm density to prevent wearing out or crushing
- Available in 4, 6 and 8mm rises

| | Code | | | Size Shoe Size | |
|-------------|-------------|-------------|---|-----------------|------|
| 4mm | 6mm | Size 8mm | | Size Silve Size | DVA: |
| FBS3016024S | FBS3016026S | FBS3016028S | S | Up to 37 | AR04 |
| FBS3016024M | FBS3016026M | FBS3016028M | M | 38 - 43 | AR04 |
| FBS3016024L | FBS3016026L | FBS3016028L | L | 44+ | AR04 |



Darco Peg-Assist

- · Designed for use in most DARCO square-toe shoes, other than the DARCO Med-Surg Shoe (for which DA-PTQ Insoles should be used)
- Removable pegs allow for localised off loading of wounds and ulcerations of the foot
- Pressure reduced by as much as 60% while allowing patient to remain ambulatory
- Poron® cover eliminates the incidence of ring oedema and edge abrasion

| Code | Size | Shoe | DVA: | |
|-------------|------|-----------|-----------|------|
| | | Female | Male | DVA: |
| FBS301699XS | XS | Up to 5 | - | AR04 |
| FBS301699S | S | 5.5 - 8 | - | AR04 |
| FBS301699M | M | 8.5 - 10 | 6 - 8 | AR04 |
| FBS301699L | L | 10.5 - 13 | 8.5 - 11 | AR04 |
| FBS301699XL | XL | 13.5+ | 11.5 - 14 | AR04 |



Rehband Adjustable Night Splint

- Patented Formfit 3D molded technology allows for variable padding thickness in strategic locations for superior comfort and fit
- Soft lycra lined foam interior is plush and breathable, eliminating heat build-up
- Three lycra lined foam padded straps with user friendly centre release buckles ensure immobilisation while maintaining comfort across the flexor tendons
- Adjustable dorsi-flexion straps provide a gentle stretching of the plantar fascia and achilles tendon.

| Code | Size | Female Shoe Size | Male Shoe Size | DVA: |
|------------|------|------------------|----------------|------|
| FBS301607 | S | Up to 7.5 | Up to 7 | AR04 |
| FBS301607M | M | 8 - 10.5 | 7.5 - 10 | AR04 |
| FBS301607L | L | 11+ | 10.5+ | AR04 |

Ankle Foot Orthoses / Toe Wiper











Rehband Hallux Valgus 510 Night Strap

- A light and flexible strap to unload the joints of the big toe
- Especially designed for nightly or post-operative usage
- Padded with skin-caring and breathable material
- Provides pain relief and guides the big toe

| Cod | le | Size | DVA: |
|------------|-----------|----------------|------|
| Left Right | | Size | DVA. |
| FBS301608L | FBS301608 | Universal Size | AR04 |



Hallux Valgus 509 Day Strap

- A light and flexible strap to unload the joints of the big toe
- Padded for optimal comfort
- Can be used with or without shoes
- Gives pain relief and may be a condition for a normalised function of the muscles, guides the big toe

| Code | | Size | DVA: |
|------------|-----------|----------------|------|
| Left Right | | Size | DVA. |
| FBS301609L | FBS301609 | Universal Size | AR04 |



Drop Foot Ankle Foot Orthosis

- Designed to hold the foot correctly
- Supports the foot to stop dragging
- Fits easily into shoe, heat mouldable, easily trimmed

| Size | Shoo Sizo | Shoe Size Code | | DVA: |
|------|------------|----------------|--------------|------|
| Size | 3110e 312e | Left | Right | DVA. |
| S | Up to 38 | FBS301589LS | FBS301589 | AR19 |
| M | 39 - 43 | FBS301589LM | FBS301589RM | AR19 |
| L | 44 - 46 | FBS301589LL | FBS301589RL | AR19 |
| XL | 47+ | FBS301589LXL | FBS301589RXL | AR19 |



Drop Foot Ankle Foot Orthosis Open

- Limits plantar flexion, common peroneal nerve injuries
- Full length sole
- · Can be adjusted for a custom fit by heating, cutting and grinding

| Size | Shoe Size | Code | | |
|------|------------|--------------|-------------|------|
| Size | Silve Size | Left | Right | DVA: |
| XS | 4 - 5 | FBS301588LXS | FBS301588 | AR19 |
| S | 5 - 7 | FBS301588LS | FBS301588RS | AR19 |
| M | 7 - 9 | FBS301588LM | FBS301588RM | AR19 |
| L | 9 - 12 | FBS301588LL | FBS301588RL | AR19 |



Toe Washer

- For washing and drying between the toes
- Long stem with plastic built up handle and a wire frame covered in a white towelling pad
- Designed to assist people with limited flexibility or access to their feet
- Includes 2 towelling pads

| Code | Description | DVA: |
|-----------|-------------|------|
| BTS113000 | Toe Washer | AU10 |

Donning / Doffing Aids







EzyAs Compression Stocking Applicator

- Unique plastic gutter to assist in pulling on compression stockings
- Stocking is fed onto gutter, held in place by side notches, slipped over foot, then pulled into place

| Code | Description | Width | Length | Calf/Arm Circ. | Warranty | DVA: |
|-----------|--------------|--------|--------|----------------|----------|------|
| DLD269001 | Small - Red | 140 mm | 140 mm | 280-350 mm | 1 Year | AU13 |
| DLD269002 | Med -Yellow | 140 mm | 140 mm | 320-400 mm | 1 Year | AU13 |
| DLD269003 | Large - Blue | 140 mm | 175 mm | 380-470 mm | 1 Year | AU13 |
| DLD269004 | Handle | 210 mm | 270 mm | - | 1 Year | AU13 |



Stocking Aid - Dorking

- A simple, inexpensive device to assist with putting on a stocking or sock especially if bending and reaching to the foot is difficult
- It consists of a strong, shaped, flexible plastic gutter with long cotton tapes attached to the top

| Code | Description | Warranty | DVA: |
|-----------|------------------------|----------|------|
| DLD269260 | Stocking Aid - Dorking | 1 Year | AU13 |



Terry Cloth Sock Aid

- Lightweight and comfortable heel contour assists changing into and out of socks
- Slits allow flexibility around the heel for best fit and two pulling tabs with loop handles permit one-handed application

| Code | Description | Width | Gutter Length | DVA: |
|-----------|----------------------|--------|---------------|------|
| DLD269090 | Terry Cloth Sock Aid | 125 mm | 230 mm | AU13 |



Stocking Aid

- Lightweight plastic gutter assists changing into and out of socks
- Sock is fed onto gutter, slipped over foot, then pulled into place

| Code | Description | Width | Gutter Length | DVA: |
|-----------|--------------|--------|---------------|------|
| DLD269275 | Stocking Aid | 130 mm | 230 mm | AU13 |



Compression Stocking Aid

- · Applies support hosiery quickly and evenly
- Grooves catch the compression stocking and prevent it from popping off
- Works like a traditional sock aid just pull the stocking over the shell, insert the foot and pull the cord to bring the hose up onto the leg
- Suitable for Class 1 and 2 stockings

| Code | Description | Warranty | DVA: |
|-----------|----------------------------------|----------|------|
| DLD269110 | Compression Stocking Aid - Small | 1 Year | AU13 |
| DLD269120 | Compression Stocking Aid - Large | 1 Year | AU13 |







Compression Stocking Aid

- EUREKA-ON Plastic Fabric Slide is a compression stocking applicator for the leg or arm
- Consists of a triangular shaped durable paper-like slippery material with a central extended pull tab and eyelet holes on the sides for threading a thin tubular plastic lacing pin

| Code | Description | Warranty | DVA: |
|-----------|-----------------------------------|----------|------|
| DLD269010 | Compression Stocking Aid - Eureka | 1 Year | AU13 |



Notch Cord Sock Aid

- Designed to assist people who have difficulty reaching their feet or the use of only one hand to put on socks or stockings
- Featuring a unique notch cut out at the lower end which hooks the sock when first putting it on the aid, making it less likely to accidentally slip off
- The attached tapes / cord are used to pull the sock aid over the feet

| Code | Description | Warranty | DVA: |
|-----------|-----------------------------|----------|------|
| DLD269100 | Sock Aid - Notch Cord Style | 1Year | AU13 |



Sock & Stocking Aid

- Designed to help those with restricted mobility, enabling the user to successfully apply on socks
- Also very effective used with stockings It consists of a plastic collar or gutter for holding the garment at one end, and a plastic shoehorn with a notch for easy sock removal at the other
- A coated steel tube connects the two

| Code | Description | Warranty | DVA: |
|-----------|--------------------------------|----------|------|
| DLD269290 | Sock & Stocking Aid - Brevetti | 1 Year | AU13 |



Tights Aid

- The Tights Aid consists of two connected plastic contoured shapes that taper at the end, with white tapes
- The foot portion of the tights are pulled onto the shapes whilst the user is seated

| Code | Description | Warranty | DVA: |
|-----------|-------------|----------|------|
| DLD269280 | Tights Aid | 1 Year | AU13 |



Easy-Slide - Applicator For Open-Toe Stockings

- · Application aid for open-toe stockings
- Also suitable for hospital stockings with inspection hole
- Slide stocking over, and then remove through open toe

| Code | Size | Shoe Size | Suit Foot | DVA: |
|-----------|---------|-----------|------------------------------|------|
| GMS340951 | Small | 34-37 | Child / Women's Small | AU13 |
| GMS340952 | Medium | 38-41 | Women's Medium / Men's Small | AU13 |
| GMS340953 | Large | 42-45 | Woman's Large / Men's Medium | AU13 |
| GMS340954 | X-large | 46+ | Men's Large | AU13 |

Donning / Doffing Aids







JOBST Stocking Donner

- Application aid for closed-toe garments
- Slide stockings over, then remove from top of stocking
- Suggest to use with rubber gloves

| Code | Size | Shoe Size | DVA: |
|-----------|-------------|-----------|------|
| PTS641090 | Small / Med | Up to 38 | AU13 |



Easy 'Sim - Slide' On / Off Aid (Open Toe)

- The ultimate donning and doffing aid for open toe compression stockings and tights
- Reduces the physical strain for both the caregiver and user
- Available in 4 sizes

| Code | Size | Shoe Size | DVA: |
|-----------|---------|-----------|------|
| PTS642508 | Small | 34 - 37 | AU13 |
| PTS642509 | Medium | 38 - 41 | AU13 |
| PTS642510 | Large | 42 - 45 | AU13 |
| PTS642511 | X Large | 46 | AU13 |



Easy Slide 'Magnide' On / Off Aid (Closed Toe)

- The easiest donning and doffing aid for closed toe compression stockings and tights
- Reduces the physical strain for both the caregiver and user

| Code | Size | Shoe Size | DVA: |
|-----------|---------|----------------|------|
| PTS642514 | 40 - 45 | Medium / Large | AU13 |



Easy Slide 'Magnide' Applicator (Closed Toe)

- The ultimate donning aid for closed toe compression stockings and tights
- Reduces the physical strain for both the caregiver and user

| Code | Size | Shoe Size | DVA: |
|-----------|------|-----------------|------|
| PTS642512 | 40 | Small / Medium | AU13 |
| PTS642513 | 41 | Large / X Large | AU13 |



Medi-Butler Stocking Donner

• The Medi-Butler is used for applying compression knee highs and thigh highs

| Code | Size | DVA: |
|-----------|---------|------|
| PTS642542 | Regular | AU13 |

Donning / Doffing Aids / Laces









Medi-Butler Off - Doffing Aid

- This long handled device assists with the removal of compression stockings, socks and normal pantyhose
- It can also double as a shoehorn

| Code | Size | DVA: |
|-----------|-----------|------|
| PTS642543 | Universal | AU13 |



JOBST Stocking Donner

- Application aid for closed-toe garments
- Slide stockings over, then remove from top of stocking
- Suggest to use with rubber gloves

| Code | Size | Shoe Size | DVA: |
|-----------|-------------|-----------|------|
| PTS641090 | Small / Med | Up to 38 | AU13 |



Sigvaris Doff N' Donner

- Soft device that allows you to quickly guide a compression garment up and down your limb
- The Cone acts as a loading device to make using the Sigvaris Doff N' Donner even easier
- Sold seperately or as a set

| Code | Description | DVA: |
|-----------|------------------------------|------|
| PTS642660 | Doff N' Donner Device | AU13 |
| PTS642661 | Doff N' Donner Cone (Black) | AU13 |
| PTS642662 | Doff N' Donner Device & Cone | AU13 |



Jobst Donning Gloves

- Jobst gloves help create friction, provide grip and help to evenly distribute the fabric on the leg
- This product contains natural rubber Latex

| Code | Description | DVA: |
|-----------|-------------------------------|------|
| FBS303490 | Jobst Donning Gloves - Small | AU13 |
| FBS303491 | Jobst Donning Gloves - Medium | AU13 |



Spring Laces

- Elasticised and spring coiled shoe laces convert lace-up shoes into slip-on shoes
- Laces fasten by wrapping the ends together so do not need to be tied into a bow

| Code | Description | Warranty | DVA: |
|-----------|-----------------------------|----------|------|
| DLD269140 | Spring Laces - Black 600 mm | 6 Months | AU14 |
| DLD269150 | Spring Laces - White 600 mm | 6 Months | AU14 |
| DLD269151 | Spring Laces - Brown 600 mm | 6 Months | AU14 |
| DLD269153 | Spring Laces - Blue 600 mm | 6 Months | AU14 |




eRAP: DVA Prescriber Portal

eRAP is Aidacare's mobile platform for Health Providers to submit RAP direct order forms. Prescribers can access our eRAP via the web portal or via our iPad app.



For more information, visit erap.aidacare.com.au or call 1300 888 052

Other Resources







1300 133 120

